
EN 4	OASIS® One Action Stent Introduction System Instructions for Use
CS 6	Jednofázový zaváděcí systém stentu OASIS® Návod k použití
DA 8	OASIS® stentindføringsystem med én funktion Brugsanvisning
NL 11	One action stentintroductiesysteem (OASIS®) Gebruiksaanwijzing
FR 13	Système d'introduction d'endoprothèse en un temps OASIS® Mode d'emploi
DE 16	OASIS® One Action Stent-Einführsystem Gebrauchsanweisung
EL 19	Σύστημα εισαγωγής ενδοπρόσθεσης μίας ενέργειας OASIS® Οδηγίες χρήσης
HU 22	OASIS® egyműveletű sztent felvezetőrendszer Használati utasítás
IT 24	Sistema di introduzione di stent in un'unica azione OASIS® Istruzioni per l'uso
NO 27	OASIS® One Action stentinnføringsssystem Bruksanvisning
PL 30	OASIS® Jednostopniowy System Wprowadzania Stentu Instrukcja użycia
PT 32	Sistema introdutor de stent com um só movimento OASIS® Instruções de utilização
ES 35	Sistema de introducción de stents en un paso OASIS® Instrucciones de uso
SV 38	OASIS® One Action stentinföringsssystem Bruksanvisning

Removable Stylet

Snímatelný stilet
 Aftagelig stilet
 Verwijderbaar stilet
 Stylet amovible
 Entfernbarer Mandrin
 Αφαιρούμενος στειλειός
 Ektávolítható mandrin
 Mandrino rimovibile
 Stilet som kan fjernes
 Usualny mandryn
 Estilite amovível
 Estilite extraíble
 Borttagbar stilet

Proximal Wire Port

Proximální vodící port
 Proksimal kateterlederport
 Proximale voerdraadpoort
 Orifice pour guide proximal
 Proximaler Drahtport
 Εγγύς θύρα σύρματος
 Proximális drótnyílás
 Ingresso prossimale della guida
 Proksimal vaierport
 Proksymalny port dla prowadnika
 Orificio proximal do fio
 Acceso proximal de la guía
 Proximal ledarport

IDE Port

Port IDE
 IDE-port
 IDE-poort
 Orifice IDE
 IDE-Port
 Θύρα IDE
 IDE-port
 Ingresso IDE
 IDE-port
 Port IDE
 Orificio IDE
 Acceso de IDE
 IDE-port

Pushing Catheter


Tlačný katetr
 Skubbekateter
 Pushing-katheter
 Cathéter pousse-prothèse
 Schiebekatheter
 Καθετήρας ώθησης
 Tolókatéter
 Catetere di spinta
 Skyvekaterer
 Cewnik popychający
 Cateter posicionador
 Catéter empujador
 Tryckkateter

Guiding Catheter

Vodící katetr
 Styrekateter
 Geleidekatheter
 Cathéter guide
 Führungskatheter
 Οδηγός καθετήρας
 Vezetőkatéter
 Catetere guida
 Ledekaterer
 Cewnik prowadzący
 Cateter guía
 Catéter guía
 Styrkateter

Radiopaque Marker

Rentgenkontrastní značka
 Røntgenfast markør
 Radiopake markering
 Marqueur radio-opaque
 Röntgengedichte Markierung
 Ακτινοσκοπικός δείκτης
 Sugárfogó marker
 Marker radiopaco
 Radioopak markør
 Znacznik cieniодajny
 Marcador radiopaco
 Marcador radiopaco
 Röntgentät markering

Proximal Wire Port

Proximální vodič port
 Proximal kateterlederport
 Proximale voerdraadpoort
 Orifice pour guide proximal
 Proximaler Drahtport
 Εγγύς θύρα σύρματος
 Proximális drótnyílás
 Ingresso prossimale della guida
 Proksimal vaierport
 Proksymalny port dla prowadnika
 Orificio proximal do fio
 Acceso proximal de la guía
 Proximal ledarport

Radiopaque Bands

Rentgenokonstránní proužky
 Røntgenfaste markeringer
 Radiopake banden
 Bandes radio-opaques
 Röntgendichte Markierungen
 Ακτινοσκοπικές ταινίες
 Sugárfogó sávok
 Bande radiopache
 Radioopake bånd
 Opaski cieniodajne
 Bandas radiopacas
 Bandas radiopacas
 Röntgentäta band

Positioning Sleeve

Zaváděcí objímka
 Positioneringsshylster
 Positioneringsshuls
 Gaine de positionnement
 Positionierhülse
 Χιτώνιο τοποθέτησης
 Pozicionáló hüvely
 Manicotto di posizionamento
 Posisjoneringshylse
 Rękaw umiejscawiający
 Manga de posicionamento
 Manguito posicionador
 Placeringshylsa

Stent

Stent
 Stent
 Stent
 Endoprothèse
 Stent
 Ενδοπρόσθεση
 Szent
 Stent
 Stent
 Stent
 Stent
 Stent
 Stent
 Stent
 Stent

Introducer

Zavaděč
 Indfører
 Introducer
 Introducteur
 Einführungsinstrument
 Εισαγωγέας
 Felvezető
 Introduttore
 Innføringsenhet
 Introduktor
 Introduitor
 Introduitor
 Införare

INTENDED USE

This device (with preloaded stent, if applicable) is intended for endoscopic biliary stent placement to drain obstructed bile ducts.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications associated with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Wire guide diameter and inner lumen of wire-guided device must be compatible.

Not compatible with pigtail stents.

Do not use this device if stent cannot advance through strictured area.

Preloaded stents should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

SYSTEM PREPARATION

1. If applicable, pre-load desired stent and positioning sleeve onto tip of introducer.

INSTRUCTIONS FOR USE

Illustrations

I. If using Intra Ductal Exchange (IDE) port and a pre-positioned short wire. (See fig. 1)

1. Unlock short wire from Wire Guide Locking Device and advance introducer with preloaded stent over pre-positioned wire guide ensuring that wire guide exits guiding catheter at IDE port.
2. Introduce device into endoscope accessory channel and relock wire guide.
3. With elevator open, advance device in short increments until endoscopically visualized exiting scope. **Note:** As front flap of stent enters accessory channel, keep positioning sleeve over back flap until completely in accessory channel. Slide positioning sleeve over pushing catheter, keeping it clear of accessory channel.
4. Advance device into appropriate duct.
5. Fluoroscopically monitor radiopaque band(s) on guiding catheter. When a sufficient length of guiding catheter is above majority of stricture, disconnect Luer lock fitting on introducer. **Important:** A sufficient portion of guiding catheter must be above stricture prior to stent positioning.
6. Maintain position of guiding catheter and advance stent into desired position using pushing catheter.
7. Fluoroscopically and endoscopically confirm desired stent placement.
8. Unlock wire guide from Wire Guide Locking Device.
9. Retract wire guide until it exits guiding catheter at IDE port.

Note: Fluoroscopically visualize radiopaque band at IDE port. When radiopaque distal tip of wire guide passes band, a disengagement from wire guide lumen will occur.

10. Advance disengaged wire guide to maintain ductal access.
11. Relock wire guide into Wire Guide Locking Device.
12. Retract guiding catheter into endoscope while maintaining position of pushing catheter.
13. Confirm position of stent once guiding catheter is removed.
14. Remove pushing catheter.

II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide. (See fig. 2)

Note: For best results, wire guide should be kept wet.

1. Remove stylet, if applicable.
2. Advance introducer with preloaded stent over pre-positioned wire guide and into accessory channel of endoscope until proximal end of wire guide exits PWP. **Note:** If applicable, ensure that wire guide does not exit at IDE port.

REFER TO STEPS 3-7 IN "SECTION I", THEN RESUME WITH STEP 3 BELOW:

3. Retract and remove both guiding catheter and wire guide into endoscope while maintaining position of pushing catheter.
4. Confirm position of stent once guiding catheter and wire guide are removed.
5. Remove pushing catheter.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2019
COOK IRELAND LTD.